



Abri-Soft

Super-absorbent and soft underpads
made in Europe



Disposable fluff underpads

Eco-labeled underpads manufactured in Europe

Abri-Soft is a range of soft and absorbent underpads that are manufactured at our production facility in France using state-of-the-art, energy-optimized machinery.

Our strong focus on sourcing the best raw materials from Europe also means that Abri-Soft has obtained several eco-labels from third parties including the Nordic Swan Ecolabel. These labels serve as your guarantee that our products are safe to use on even the most sensitive skin.

Through routine quality control and in-house quality testing, we strive to ensure that Abri-Soft continues to be one of the best underpads on the market.



The best possible protection

All products in our bed protection range are composed of different layers and materials. Combined, they provide a high-absorbing underpad that feels comfortable on the user's skin.

Abri-Soft technology:

- Sealed barriers and a waterproof PE back sheet for optimal leakage protection
- A soft, non-woven top layer and a crackle-free material ensures comfort for the user
- A fast-absorbing core made of long-fiber fluff (and SAP in the Superdry range) ensures a dry-to-the-touch surface
- An even distribution of fluff (and SAP) ensures a thin and discreet product
- Smart packaging makes it easy to grab one underpad at a time to prevent cross-contamination

An absorbing underpad for every need:

- 7 ranges of absorbencies
- Wide variety of different sizes

By choosing the right product for your needs, you minimize the number of products used and help reduce waste.



Ultra Light



Light



Basic



Superdry



Classic



Premium



Excellent



More than incontinence:

A product for every need

While often known for their use as extra security for people living with incontinence, the great variety of sizes and absorbencies of the Abri-Soft range make it ideal for several areas of use.

Ideal protection anywhere you need to absorb liquids:



Extra protection for incontinence



Bed bathing and hygiene



Wound care



Surgery



Birth



Examination and clinic



Blood tests

Certified by third parties

FSC™ (Forest Stewardship Council) is an international non-governmental organization dedicated to the promotion of responsible forest management worldwide. FSC™-certified products are made from raw materials from well-managed forests and responsible sources. The label's traceability and integrity are upheld through stringent control measures, and by choosing FSC™-certified products, you help take care of forests and the people and wildlife that call them home. For more information, go to www.fsc.org.

The Swan label, also known as the Nordic Swan Ecolabel, is an official Scandinavian environmental label. Products with the Nordic Swan Ecolabel meet stringent environmental requirements with regard to raw material production, use, documentation, and safety.

The Dermatologically Tested symbol guarantees that products have been tested for potential skin irritations by an independent third-party laboratory.

About ABENA

ABENA, founded in 1953, is a Danish, family-owned manufacturer and wholesaler of health care solutions and related products headquartered in Southern Denmark.

The ABENA Group, which is represented by subsidiaries in a long list of countries, supplies more than 36,000 products and operates in more than 90 countries across the globes.

ABENA-FRANTEX is based in Nogent sur Oise, north of Paris. ABENA-FRANTEX has been a part of the ABENA group since 1992, but has been producing absorbent underpads for the ABENA group since 1988.

BECAUSE WE CARE